

Humidifier Adapter Tubing, 12"

HCPCS:



- The Humidifier Adapter Tubing allows for patients to attach their humidifier bottle to their oxygen concentrator for humidified oxygen therapy.
- Its three-channel crush resistant design prevents blockages and provides superior safety.
- Attaches humidifier bottle to oxygen concentrator for more comfortable, humidified oxygen therapy
- Three-channel safety oxygen supply tube with one molded trumpet connector end and one male thread grip end
- Durable clear tubing that is crush and kink resistant
- 12" long for close connections
- Compatible with most humidifier bottles and oxygen concentrators
- Quantity: 25/case

PRODUCTS

Description

9996-1	Humidifier Adapter Tubing, 12 25EA/CV
--------	---------------------------------------