

M6 Cylinder Rack, 12 Cylinder Capacity

HCPCS:

- The M6 Cylinder Rack ensures durability and strength in storing your oxygen cylinders with an all-welded frame construction.
- It also features individual spacing, positive cylinder placement, and weight distribution to promote safe storage.
- Features an all-welded frame construction with crossbars to ensure durability and strength
- Promotes safe storage with individual spacing, positive cylinder placement, and weight distribution
- Equipped with bolt down bars and feet
- Chrome finish
- Fits up to 12 M6 cylinders
- Sold individually

PRODUCTS

Description

CR-M6PC12

M6 OXYGEN CYLINDER RACK-12 CYL