

D or E Cylinder Cart, Silver Vein

HCPCS:



- The Roscoe Adjustable Oxygen Carts make it easy to transport your oxygen tank anywhere.
- They are constructed of durable, high-strength steel and can accommodate either "D" or "E" size cylinders.
- Made to last with durable, high-strength steel
- Handle is height-adjustable to accommodate most individuals
- Large, easy-roll wheels allow for safe maneuvering
- Includes gurney/mask hook, nonslip hand grip and one-piece thumb screw
- Silver vein finish; also available in chrome plating
- Can accommodate D or E cylinders
- Quantity: 4/case
- Limited lifetime warranty
- HCPCS Code: E1355*

PRODUCTS

Description

E-CARTPC

OXYGEN D/E CYLINDER CART 4EA/CV