

Chad Oxymizer Pendant - Disposable Oxygen Conserver/Reservoir

HCPCS:



- Oxymizers are Disposable Oxygen Conserver that are available as a Pendant Oxymizer or a Mustache Oxymizer.
- The Pendant Oxymizer can be used with a continuous oxygen flow of up to 15 LPM.
- This Pendant Oxymizer features a reservoir that can be worn under clothing.
- The Mustache Oxymizer can be used with a continuous oxygen flow of 20 LPM.
- Both Oxymizers provide continuous high-flow oxygen therapy for homecare, hospice, clinic, hospital, or long-term care.
- The Oxymizers are compatible with oxygen concentrators, compressed oxygen cylinders, and liquid oxygen.
- Oxymizers can be used with pursed lip breathing.
- Oxymizers may be used during exercise, resting, or during sleep and offer patients a comfortable alternative to using an oxygen mask.
- Features & Benefits
- Disposable Nasal Cannula Oxygen Conserver for Oxygen Savings
- Pendant Oxymizer features a reservoir that can be worn under clothing

PRODUCTS	Description
P-224EA	PENDANT STYLE OXYMIZER 1EA/