

Portable Lightweight Transport Wheelchair with Detachable Footrests, Midnight Blue

HCPCS:



FIGURE A

- Lightweight transport wheelchairs by Hugo make it easier than ever for you to help your loved ones get around in comfort.
- This portable transport wheelchair is much lighter than standard transport chairs at only 16 lbs (without accessories) and ultra-compact when folded.
- Easy-touch levers allow you to fold the back of the chair for more compact storage. The adjustable, detachable footrests make this mobility aid great for traveling and easily transportable in a car trunk, bus, plane, etc.
- Ergonomically angled handlebars increase comfort, offering exceptional shoulder room for both the caregiver and 35% more shoulder room for the user, making it easier to push.
- Large 8" wheels provide increased comfort when overcoming obstacles. The Hugo Transport Wheelchair has a fully padded seat and backrest, as well as a convenient storage pocket.
- A seat belt is included as a standard safety feature while the dual locking rear wheels provide for secure transfers.
- Tool-free assembly and durable yet lightweight aluminum frame that supports up to 250 lbs make Hugo Transport Wheelchairs a great mobility aid choice.

PRODUCTS

Description

700-869

Hugo TC, KD Frame, Blue, 1 EA,2CV,Ret