

Pure Expression Carry Bag

HCPCS:



- Carry bag for RTLBP2000

PRODUCTS	Description
BP001	Breast Pump Carry Bag,1/ea
BP0011	PowerAdapter,6V,RTLBP2000/RTLBP2200, 1e
BP0012	ExternalBatteryPack,6V,RTLBP2000/BP2200
BP0013	PumpKit w/o Insert,31MM,BP1000/2000/220
BP0014	Silicone Insert, 27MM, BP1000/2000/2200
BP0015	FeedingNippleSet, BP1000/2000/2200,1/pk
BP0016	Base for Bottles, BP1000/2000/2200, 1/e
BP002	Insulated Cooler Bag, 1/ea
BP005	PumpKit w/o Insert,27MM,BP1000/2000/220
BP006	Silicone Insert, 24MM, BP1000/2000/2200
BP007	ValveSetW/Membrane,BP1000/2000/2200,1/p
BP008	Membranes ,BP1000/2000/2200, 4/bag
BP009	Tubing Set, BP1000/2000/2200, 1/pr