



## 800cc Hi-Flow Suction Canister with Aerostat Filters

HCPCS:

- These 800cc Hi-Flow Suction Canisters are a high-quality way to meet fluid collection and disposal needs.
- They provide a high flow rate, so suctioning is quick and efficient during procedures.
- For small-volume use
- Includes self-sealing lid with leak-free seal
- Features pour spout for emptying without removal of lid, vacuum and horizontal patient port for direct connection of tubing
- Features Aerostat bacterial filter to protect against contamination
- Built-in critical measure for more accurate measuring
- Mechanical shutoff valve prevents overflow and contamination
- Universal ports and size for use with most suction machines
- Disposable, single use
- Dimensions: 6" H x 9" L x 4" W

### PRODUCTS

### Description

1160B-1

BEMIS 8002055 800ML SUCTION 1EA