

Bulb-Tip Yankauer with 6' Tubing

HCPCS:

- The Bulb-Tip Yankauer allows for precise suctioning without damaging the surrounding tissue.
- It features a slip-resistant handle for high-quality performance and comes equipped with 6-foot tubing.
- Provides precise and effective suctioning during medical procedures
- Comes equipped with 6-foot tubing
- Does not include control vent
- Features slip-resistant handles, smooth inner surface and uniform tip-to-base internal diameter for fast evacuation
- Clear, transparent construction for easy visualization of fluids
- Ribbed five-in-one connector; fits various sizes of suction tubing
- Sterile, for single patient use
- Quantity: 20/case

PRODUCTS

Description

DYND50135

Bulb-Tip Yankauer w/6 ft. 1CV/