

Suction Setup Kit

HCPCS:

- The Suction Setup Kit allows users to conveniently replace all necessary parts in their suction machine at once.
- It includes a yankauer with bulb tip and control vent, replacement elbow, bacteria filter, 800ml reusable collection jar, 6' tubing and 16" vinyl tubing.
- Each item in the kit can also be purchased separately if desired:
- Bulb-tip Yankauer - #55K86V
- Elbow - #ROS-610E
- Bacteria Filter - #ILB-38B
- Collection Jar - #1160B-1
- 6' Tubing - #55N56A
- 12 Tubing #SUC12-1
- To keep your portable suction machine performing optimally, replace its parts regularly and easily with the Suction Setup Kit.

PRODUCTS

Description

KIT-018D	SUCTION SET UP KIT- DEVILBISS 1EA/
----------	------------------------------------