

REPLACEMENT HAND NUT FOR 1EA/

HCPCS:

PRODUCTS	Description
1160B-1	BEMIS 8002055 800ML SUCTION 1EA
55K86	YANKAUER SUCTION HANDLE W/O VENT 50EA/CV
55K86V	YANKAUER SUCTION HANDLE W/VENT 50EA/CV
ASP-KIT	ASP-ROSP SUCTION KIT 12PKG/CV
DYND50135	Bulb-Tip Yankauer w/6 ft. 1CV/
KIT-018D	SUCTION SET UP KIT- DEVILBISS 1EA/
KIT-019	SUCTION SETUP KIT 1EA/
KIT-020	SUCTION SETUP KIT 1EA/
KIT-50004	Suction Kit for 50004 & 50006, 1PKG/
KIT-7305D	Suction kit for 7305D & 7305P, 1PKG/
TK-SUC4	Silicone Tubing for Suction 4 1PKG/