

## SOFT PAD RAISED TOILET SEAT

HCPCS:



| PRODUCTS      | Description                              |
|---------------|------------------------------------------|
| 90077         | Bracket for Toilet Safety 1EA/           |
| FGB302C0 0000 | RAISED TOILET SEAT 1EA/                  |
| FGB30300 0000 | E-Z LOCK RAISED TOILT SEAT W/ARMS 2EA/CV |
| FGB30400 0000 | E-Z LOCK RAISED TOILT SEAT W/ARMS 2EA/CV |
| FGB30500 0000 | E-Z LOCK RAISED TOILET SEAT 2EA/CV       |
| FGB30600 0000 | TOILET SEAT ELEVATOR ELONGATED 4EA/CV    |
| FGB30700 0000 | TOILET SEAT ELEVATOR ROUND 4EA/CV        |
| FGB31000 0000 | RAISED TOILET SEAT - RETAIL 2EA/CV       |
| FGB311C0 0000 | E-Z LOCK RAISED TOILET SEAT W/ARMS       |
| FGB312C0 0000 | E-Z LOCK RAISED TOILET SEAT 1EA/         |
| FGB31300 0000 | SAFE LOCK RAISED TOILET SEAT 3EA/CV      |
| FGB313CA      | SAFE LOCK RAISED TOILET SEAT 3EA/CV      |
| FGB31600 0000 | TOILET SEAT ELEVATR W/HNDLS ELONG 2EA/CV |
| FGB31700 0000 | TOILET SEAT ELEVATR W/HNDL STNDRD 2EA/CV |
| FGB31877 0000 | RAISED TOILET SEAT W/ ARMS 3EA/CV        |
| FGB32000 0000 | LOCKING RAISED TOILET SEAT 2EA/CV        |
| FGB32100 0000 | HINGED TOILET SEAT RISER ELONG 3EA/CV    |
| FGB32200 0000 | HINGED TOILET SEAT RISER STANDARD 3EA/CV |
| FGB82100 0000 | SOFT PAD RAISED TOILET SEAT              |