

ProBasics Aluminum Two-Button Release Folding Walker Without Wheels, 4/CS

HCPCS:

- ProBasics Two-Button Release Folding Walkers provide maximum strength in a lightweight aluminum design, all backed by a limited lifetime warranty.
- Ample height adjustment options coupled with a two-button folding mechanism make ProBasics Folding Walkers easy to use and accommodating to the most demanding patient needs.
- Available with or without 5" wheels and in junior or adult sizes.
- Support patient weights up to 300 pounds.
- Attractive, lightweight aluminum frame
- Adjustable height settings
- Easy-to-use push button folding mechanism provides audible click when locked
- Front cross brace near top of walker allows for full stride
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Vinyl contoured hand grip for added comfort and long lasting wear
- Rear glide cap allows for walker to slide easily and smoothly over most surfaces (come standard on models with 5" wheels)

PRODUCTS

Description

WKAAN2B

PROBASICS ALUM ADULT WALKER, 2 4EA/CV