

## Electrode - Poly - 2x3.5 -

HCPCS:

- The White Cloth 2" x 3.5" Electrodes are made with the highest clinical standards and are for use with TENS and EMS units.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- They are self-adhesive and easy to use, and with premium white cloth and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort on a variety of application sites
- Durable and comfortable premium white cloth backing
- Durable and reliable full-coverage wire connector
- 4 electrodes/pack; comes sealed in tamper-proof poly bag
- Size: 2" x 3.5"
- Can be reused up to 15 times per pack
- Purchase the White Cloth 2" x 3.5" Electrodes for effective yet economic electrotherapy.

### PRODUCTS

### Description

E1P2035WC2

Electrode - Poly - 2x3.5 - 320EA/CV