

Electrode - Foil - 1.25 -

HCPCS:

- The 1.25" Round Tan Cloth Electrodes are made with the highest clinical standards and are for use with TENS and EMS units.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- They are self-adhesive and easy to use, and with premium tan cloth and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort
- Durable and comfortable premium tan cloth backing
- Durable and reliable full-coverage wire connector
- 4 electrodes/pack; comes sealed in tamper-proof foil bag
- 10 packages/foil bag
- Size: 1.25" round
- Can be reused up to 15 times per pack

PRODUCTS

Description

EF1200TC2	Electrode Foil 1.25" 1EA/
-----------	---------------------------
