

InTENSity Self-Adhesive Electrodes, 2" x 2" White Cloth

HCPCS:

- The InTENSity 2" x 2" square electrodes are made with the highest clinical standards and fit most TENS, EMS, IF, and other electro-stim applications.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- They are self-adhesive and easy to use, and with premium white cloth and a full-coverage wire connector, InTENSity electrodes are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort on a variety of application sites
- Durable and comfortable premium white cloth backing
- Durable and reliable full-coverage wire connector
- 4 electrodes/pack
- Can be reused up to 15 times per pack
- Purchase InTENSity Electrodes for effective yet economic electrotherapy.

PRODUCTS	Description
EP2020WC2-INTM	InTENSity Mini Pack, WC, 2 x 2 640PKG/CV