

TM DI1010 KIT 1EA/

HCPCS:

PRODUCTS	Description
77645	INTELECT TENS PURCHASED 1EA/
77711	INTELECT TENS BASIC PURCHASED 1EA/
77713	INTELECT TENS DIGITAL PURCHASE 1EA/
ACRL-2001	ACCURELIEF SINGLE CHANNEL 6EA/CV
ACRL-3001	ACCURELIEF DUAL CHANNEL 6EA/CV
DA1812	TENS AA - Digital TENS Unit 40EA/CV
DI0007	INTENSITY 7 OTC 32EA/CV
DI1010	3nTENSity 10 Digital TENS 32EA/CV
DI1212	InTENSity 12 32EA/CV
DI3302	InTENSity 5000 - Hybrid TENS 32EA/CV
DT2020	Venti TENS w/ 12 preset 20EA/CV
DT3002	TENS 3000 Analog 32EA/CV