

## INTENSITY 7

HCPCS:



- The InTENSity 7 is one of the most popular, and best valued, digital TENS units on the market and is now available without a prescription.
- This digital TENS unit is proven to be among the most durable and easy-to-use TENS units providing drug-free pain relief for sore or aching muscles of the lower back, shoulder, arm, leg due to strain from exercise and normal activities.
- The upgraded digital processor is more precise and accurate than ever before.
- This unit is great for users seeking a digital device that is also equipped with physical amplitude knobs.
- FDA cleared for over-the-counter sale.
- Providing Drug-Free Pain Relief for Sore or Aching Muscles of the Lower Back, Shoulder, Arm, Leg Due to Strain From Exercise and Normal Activities
- Modes: Burst, Normal, Modulation, SD1 and SD2
- Dual Channel for Use with Two or Four Electrodes
- Pulse Width and a Pulse Frequency Are Adjustable
- System Includes Device, Lead Wires, Four Self-Adhesive Reusable Electrodes, 9-Volt Battery, Hard Plastic Carrying Case, Quick Start Guide and Instruction Manual
- Compliance Meter, 60 Records Up to 999 Hours

| PRODUCTS | Description             |
|----------|-------------------------|
| DI0007   | INTENSITY 7 OTC 32EA/CV |