

Venti TENS Digital Pain Relief System

HCPCS:



- The Venti TENS by Roscoe Medical is a digital TENS device with three waveforms and twenty treatment modes to choose from.
- This incredibly advance, pocket sized portable device is known for its body part specific settings to help patients with easy set-up and application of each program.
- The flexibility of the treatment programs (12 presets and 8 manual modes) provides relief for both acute and chronic pain conditions.
- Dual Channel (4-lead) digital TENS
- Soft comfortable design with LCD Display
- 12 preset and 8 manual treatment modes
- Three waveforms
- Symmetrical Biphasic
- Asymmetrical Biphasic
- Monophasic
- Body part diagram for targeted pain relief

PRODUCTS

Description

DT2020

Venti TENS w/ 12 preset 20EA/CV