

InTENSity Select Combo II

HCPCS:

- The InTENSity Select Combo II is the newest standard in portable, multifunction electrotherapy offering patients four (4) modalities in one user-friendly, ergonomic design - Transcutaneous Electrical Nerve Stimulation (TENS), Neuromuscular Electrical Stimulation (NMES), Interferential Stimulation (IF), and RUSSIAN (RUSS).
- The wide-ranging capability of the InTENSity Select Combo II provides the most up-to-date treatment technology and waveforms.
- Equipped with pre-set body part specific settings as well as a customizable option, it intuitively guides the patient step-by-step to ensure the ideal therapeutic treatment throughout the healing cycle.
- This powerful device delivers up to 10,000 Hz of IF therapy with sweep frequency capabilities to prevent accommodation while providing maximum treatment benefits.
- Simply choose a preset program by body part or customize a program to meet your specific therapy needs.
- The RUSSIAN setting provides 2500Hz for a deeper and more comfortable stimulation helping to reduce muscle spasms as well as muscle strengthening.
- Contraction and relaxation times are adjustable.
- With the powerful rechargeable lithium ion battery, this portable device can be used at home, at work or on the go.

PRODUCTS	Description
DI2195	InTENSity Select Combo II 32EA/CV