

## TM DS2202 KIT 1EA/

HCPCS:

PRODUCTS	Description
77716	INTELECT NMES DIGITAL PURCH 1EA/
ACRL-5500	ULTIMATE FOOT CIRCULATOR 2EA/CV
ACRL-9001	ACCURELIEF DUAL WIRELESS 6EA/CV
ACRL-9100	ACCURELIEF WIRELESS 6EA/CV
DE5030	EMS 5.0 - Analog EMS 32EA/CV
DE7502	EMS 7500 32EA/CV
DI2195	InTENSity Select Combo II 32EA/CV
DI2717	InTENSity Twin Stim IV 32EA/CV
DI2738	InTENSity IF Combo II 32EA/CV
DI3717	InTENSity Twin Stim 32EA/CV
DI4738	InTENSity IF Combo (TENS/IF) 32EA/CV
DI8195	InTENSity Select Combo - 32EA/CV
DQEMG	EMG PRO WITH ELECTRODES 1EA/
DQEMGP	EMG PRO + PROBES 1EA/
DS2202	Twin Stim 2nd - TENS/EMS 32EA/CV
DS5002	TwinStim Plus 2nd Edition 20EA/CV
DS5402	Twin Stim Plus 3rd Edition 1EA/
DT4222	CareTec II - Digital TENS/EMS 32EA/CV
DT4444	CareTec IV - Digital 32EA/CV