

Twin Stim Plus 3rd Edition

HCPCS:

- The Twin Stim Plus 3rd Edition is a four (4) channel device that combines TENS, EMS, IF (Interferential) and Russian therapies in a single portable unit.
- With separate intensity controls for each channel and fully adjustable parameters, patients can easily customize their treatment.
- In addition, with four channels, patients have the ability to treat large areas (entire back) or multiple treatment sites at the same time.
- Other features include adjustable ramp time, large LCD screen display, two treatment timers and patient compliance meter to store operation records.
- Simple to operate, portable device for pain management
- Reduces inflammation and edema
- Great for muscle re-education
- Increases local blood circulation and relaxes muscle spasms
- Intensity level is controlled by the press of a button
- Four Stimulation Programs: TENS, EMS, IF and RUSS (Russian)
- Five modes for TENS therapy: Burst, Normal (N), Modulation, Strength-Duration 1 and 2

PRODUCTS

Description

DS5402	Twin Stim Plus 3rd Edition 1EA/
--------	---------------------------------
