

## Relief Wrap Conductive Brace with 2" x 2" Electrodes, Medium

HCPCS:

- The Relief Wrap Conductive Brace helps users receive electrode therapy in hard to reach areas without requiring the help of others.
- It can be customized to fit patient's needs based on treatment protocol, accommodating either 4 electrodes for TENS and EMS therapy or 8 electrodes for interferential therapy.
- Comes with 4-pack of attachable 2" x 2" electrodes (#EP2020RW2)
- Size medium, fits approximately 25" to 30" waists; also available in small, large, and extra large
- Easy to maintain, hand washable
- Compatible with most devices

### PRODUCTS

### Description

---

GF5020	Relief Wrap - Medium 2x2 1EA/
--------	-------------------------------

---