

Garmetrode Conductive Glove - Universal Fit

HCPCS:

- The Garmetrode Conductive Glove provides 360 degrees of even electrode treatment over the hand and wrist for the utmost effectiveness in treating injuries and disorders.
- It is compatible with all standard TENS lead wires pins and its stretchable fabric allows for a comfortable, universally snug fit.
- Effective and convenient, providing even stimulation to the entire hand and wrist
- Universal size with stretchable fabric
- Compatible with standard TENS lead wires pins (0.080 in. diameter)
- Reusable and hand washable
- Includes: Four (4) EF3000TC2, 3" round electrodes and one (1) WS3147 pig tail snap

PRODUCTS

Description

GU2856	Garmetrode Glove - Universal 1EA/
--------	-----------------------------------
