

Sliding Transfer Bench

HCPCS:

- The Sliding Transfer Bench with Swivel Seat provides an easy way to transfer in and out of the tub.
- The swivel seat turns 360° and locks every 90° allowing for easy entry and exit from the bench.
- The seat slides smoothly left and right along the frame, and comes equipped with a seat belt for added safety.
- Swivel seat turns 360° and locks every 90° allowing easy entry to and exit from the bench
- Specially-designed sliding rings lock the seat at each end of the sliding rails
- Seat glides over polished, high-strength aluminum tubes
- Molded plastic seat and back, textured finish
- Rust-proof, lightweight aluminum construction
- Adjustable height to fit user and bathtub clearance
- Comes with safety belt strap
- Sold as each

PRODUCTS	Description
77662	Sliding Transfer Bench 1EA/