

## Wand for Soundcare Plus, 1cm Head

HCPCS:

- This is the replacement wand for the SoundCare Plus Ultrasound Machine.
- It offers pulsed or continuous ultrasound therapy at 1MHz or 3MHz and includes a 1cm soundhead applicator and attached cable.
- Offers pulsed or continuous ultrasound therapy at 1MHz or 3MHz, providing optimum therapy for a range of injuries
- 1 cm soundhead applicator; also available with 5cm soundhead applicator
- Pre-calibrated, water resistant, and ergonomically designed
- Includes wand with head and cable attached
- Sold individually

| <b>PRODUCTS</b> | <b>Description</b>               |
|-----------------|----------------------------------|
| DQ9275W1        | SoundCare Plus - Wand, 1 cm 1EA/ |