

Roscoe Reusable Hot/Cold Gel Pack (5" x 10")

HCPCS:



- The Roscoe Reusable Hot/Cold Gel Pack can be heated in the microwave to help relax sore muscles and promote increased blood flow or chilled in a freezer to help reduce swelling and blood flow.
- It is made of a durable, gentle material and can be reused whenever needed.
- Can be used to treat a variety of injuries with cryotherapy or heat therapy, as it can be both chilled in a freezer or heated in a microwave
- Made of soft touch material for durability and comfort
- Stays flexible even when frozen, so pack can be molded to any body part for targeted therapy
- Reusable, can be cleaned hundreds of times
- Size: 5" x 10"
- For a versatile gel pack capable of both cryotherapy and heat therapy, purchase your Roscoe Reusable Hot/Cold Gel Pack.

PRODUCTS	Description
BG5010	Hot or Cold Gel Pack - 5x10 20EA/CV