

Roscoe Reusable Hot/Cold Gel Pack (7.5" x 11")

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- The Roscoe Reusable Hot/Cold Gel Pack can be heated in the microwave to help relax sore muscles and promote increased blood flow or chilled in a freezer to help reduce swelling and blood flow.
- It is made of a durable, gentle material and can be reused whenever needed.
- Can be used to treat a variety of injuries with cryotherapy or heat therapy, as it can be both chilled in a freezer or heated in a microwave
- Made of soft touch material for durability and comfort
- Stays flexible even when frozen, so pack can be molded to any body part for targeted therapy
- Reusable, can be cleaned hundreds of times
- Size: 7.5" x 11"
- For a versatile gel pack capable of both cryotherapy and heat therapy, purchase your Roscoe Reusable Hot/Cold Gel Pack.

PRODUCTS

Description

BG7511	Hot or Cold Gel Pack - 7.5x11 12EA/CV
--------	---------------------------------------