

## TheraMED Professional Oversized Foam-Filled Terry Cover (24" x 32")

HCPCS:

- The TheraMED Professional Oversized All-Terry Cover insulates heat to provide comfort and protect TheraMED Moist Heat Pack users from overexposure.
- It is designed with hook and loop closures to keep hot pack secure and prevent unwrapping or slipping.
- Constructed with a polyurethane foam layer between two layers of comfortable terry for optimum insulation; also available with only terry (HPC2432-TC)
- Hook and loop closures keep hot pack secure and prevent unwrapping or slipping
- Machine washable
- Oversized 24" x 32"; available in a variety of sizes
- Ideal for both high-volume clinical use or at-home therapy
- One-year warranty

### PRODUCTS

### Description

---

HPC2432-FF	Oversized Terry Cover, Foam 10EA/CV
------------	-------------------------------------

---