

AquaSense Ergonomic Adjustable Bath Seat without Backrest

HCPCS:



FIGURE A

Enjoy your bath or your shower seated on your ergonomic adjustable bath seat, by AquaSense®, a bath seat with a comfortable contoured seat! This contoured adjustable bath seat without back features a textured surface, drainage holes, seat handles, and slip resistant rubber safety tips. Adjustable legs raise or lower the seat height.

- Easily set down and access your shower spray with built in cut-outs in the seat.
- Both lightweight and durable the rust-resistant aluminum frame supports up to 300 lbs.
- Each leg includes a non-skid, non-slip tip offering security on wet surfaces.
- Simply and quickly adjust the height with the help of a push pin to achieve the ideal seated position. Seat height adjustments on this model range from 15" to 22".
- Push button, snap together design requires no tools and assembles in just minutes.
- Feel good knowing that plastics parts of this product are made using recycled materials.

PRODUCTS

Description

770-534	AQUA BATH SEAT W/O B.SUN., 1EA,2CV,RET
770-535	AqS BSeat w/o Back, White, 1EA,2CV,Ret

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Weight Capacity	300 lb