

Deluxe Bariatric Shower Chair with Cross-Frame Brace

HCPCS: E0245



- Frame has width/depth cross braces attached with aircraft-type rivets
- Durable, blow-molded plastic bench
- Aluminum frame is lightweight, durable and corrosion proof Drainage holes reduce slipping
- Euro-style clip for leg height adjustment Height adjusts in 1" increments (Figure A)
- Weight Capacity 500 lbs
- Made using recycled materials

PRODUCTS

Description

| | |
|-----------|------------------------|
| 12021KD-1 | Bath BenchHD,w/Bk,1/ea |
| 12022KD-1 | Bench,HD,w/o.Bk,1/ea |

SPECIFICATIONS

Description

| | |
|------------------------|-----------------------------|
| Warranty | Limited Lifetime |
| Seat Dimensions | 20" (W) x 12" (D) |
| Weight Capacity | 500 lb |
| Carton Shipping Weight | 32.5 lb 8.5 lb 27 lb 8.0 lb |
| Weight | 7.5 lb 6.0 lb |

ACCESSORIES

Description

| | |
|--------|--------------------------|
| 12021T | RubberTip,12021,1/ea, EA |
|--------|--------------------------|