



Shower Chair with Back and Removable Padded Arms

HCPCS:



FIGURE A

- Removable padded armrests provide comfort and added security
- Blow-molded bench and back provide comfort and strength
- Drainage holes in bench reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Leg height adjusts in 1" increments
- Made using recycled materials

PRODUCTS

Description

PRODUCTS	Description
12445KD-1	KD Bath Bench w/Back & Padded Arm,1/ea

SPECIFICATIONS

Description

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	350 lbs
Carton Shipping Weight	21 lb
Weight	9.75 lb
Width Between Arms	18.5"