

## Deluxe Handheld Shower Massager with Three Massaging Options

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Adjustable massaging head
- Chrome shower-arm mount and extra-long 84" chrome/mylar reinforced tangle-free hose
- Can be used as handheld shower spray or used in handle mount as traditional showerhead
- Attractive retail packaging
- Useful when shower, chair, bench or stool is being used

### PRODUCTS

### Description

RTL12045	Shower Massager,Deluxe,1 EA,24/CS
----------	-----------------------------------

### SPECIFICATIONS

### Description

Warranty	Limited Lifetime
----------	------------------