

ProBasics Bariatric Drop-Arm Commode,

HCPCS:



- The ProBasics Bariatric Drop-Arm Commode features a durable steel frame with an easy-release push button arm mechanism allowing for safe lateral patient transfers to and from the commode.
- Equipped with wider seat and reinforced steel frame to support patient weights up to 650 pounds, it comes complete with commode bucket and lid.
- Features a reinforced steel frame and an extra wide platform seat designed for comfort and stability
- Equipped with an easy-to-release arm mechanism allowing for safe lateral patient transfers to and from the commode
- Bucket easily slides into place
- Contoured armrests provide additional patient comfort
- Comes complete with a commode bucket and cover
- Bucket capacity: 9.5 quarts
- Color: gray frame, gray accessories
- HCPCS Code: E0168*

PRODUCTS

Description

BSBDAC

BARIATRIC DROP ARM COMMODOE 1EA/