

ProBasics Bariatric Transfer Bench, 500 lb.  
Weight Capacity

HCPCS:



- The ProBasics Bariatric Transfer Bench helps make transfers in and out of the bathtub safer and more comfortable.
- Constructed with an aluminum anodized frame, it features added cross braces for stability and a heavy duty blow-molded seat.
- It converts easily for left or right hand transfers.
- The angled legs have one inch,easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety.
- Tool-free assembly and limited lifetime warranty are standard.
- Aluminum frame is lightweight, durable and corrosion-proof
- Durable blow-molded bench and back have a textured surface and are contoured with drainage holes to reduce slipping
- Wide-based tips and suction cups add sturdiness and prevent slipping
- Legs have one-inch easy snap-button adjustments
- Added cross brace for stability
- Converts easily for left or right hand transfers

PRODUCTS

Description

|       |                                 |
|-------|---------------------------------|
| BSBTB | BARIATRIC TRANSFER BENCH 2EA/CV |
|-------|---------------------------------|