

ProBasics Padded Transfer Bench

HCPCS:



- The ProBasics Padded Transfer Bench helps make transfers in and out of the bathtub safer and more comfortable.
- Constructed with an aluminum anodized frame and padded back and bench, it features a reversible back for right or left handed transfers.
- The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety.
- Tool-free assembly and limited lifetime warranty are standard.
- Aluminum frame is lightweight, sturdy and corrosion-proof
- Vinyl seat and back are padded with 1.25" open-cell foam for added comfort
- Wide-based tips and suction cups add sturdiness and prevent slipping
- Legs have one-inch easy snap-button adjustments
- Converts easily for left or right hand transfers
- Quick tool-free assembly
- HCPCS Code: E0247*

PRODUCTS

Description

BSTBP	PADDED TRANSFER BENCH 2EA/CV
-------	------------------------------