

ProBasics Raised Toilet Seat with Lock and Arms,

HCPCS:



- ProBasics Raised Toilet Seats are a must for anyone with bending or sitting difficulties.
- Constructed of heavy duty blow-molded plastic that is durable and easy to clean.
- Adjustable front locking mechanism tightens easily with one hand and fits most toilets.
- Prevents shifting of seat when transferring.
- Padded arms for added security and comfort
- Elevates existing seat height by 4.5"
- Extra-large hygienic hole cutout
- Tool-free assembly
- Color: white
- Seat width/height: 17"/17"
- Width between arms:19.25"

PRODUCTS

Description

BSRTSLA	RAISED TOILET ST W/LOCK&ARMS 1EA/
---------	-----------------------------------