

## AquaSense 3-in-1 Contoured Raised Toilet Seat

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Whether you're post-surgery or providing your family members with independence in the washroom, this raised toilet seat will provide added support for getting safely on and off one's toilet
- 3 easy-to-assemble configurations for maximum versatility
- Use without handles for maximum space and easy access
- Use one handle for easy lateral access and a secure grip
- Use both handles for more security and support getting up and down
- Allows the user to position the armrests in order to adapt the product to his / her preference and to the different stages of recovery
- Proprietary body-contoured shape with an improved perinea scoop
- Attaches securely to most standard and elongated toilets for superior flexibility
- Raises toilet seat height 4" (10 cm) on 770-618 and 2" (5 cm) on 770-617 for better access and security
- Easy-to-clean polyethylene for fuss-free clean ups
- Corrosion resistant for increased durability
- Easy and quick set-up
- AADL# B230

### PRODUCTS

### Description

770-617	AqS 3-in-1 Raised Toilet Seat, 2
770-618	AQS 3IN1 RTS 4", WHITE, 1EA,2CV,Ret

### SPECIFICATIONS

### Description

Weight Capacity	300 lb( 136 Kg)
-----------------	-----------------