

AquaSense 3-in-1 Contoured Raised Toilet Seat

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Whether you're post-surgery or providing your family members with independence in the washroom, this raised toilet seat will provide added support for getting safely on and off one's toilet
- 3 easy-to-assemble configurations for maximum versatility
- Use without handles for maximum space and easy access
- Use one handle for easy lateral access and a secure grip
- Use both handles for more security and support getting up and down
- Allows the user to position the armrests in order to adapt the product to his / her preference and to the different stages of recovery
- Proprietary body-contoured shape with an improved perinea scoop
- Attaches securely to most standard and elongated toilets for superior flexibility
- Raises toilet seat height 4" (10 cm) on 770-618 and 2" (5 cm) on 770-617 for better access and security
- Easy-to-clean polyethylene for fuss-free clean ups
- Corrosion resistant for increased durability
- Easy and quick set-up
- AADL# B230

PRODUCTS

Description

770-617	AqS 3-in-1 Raised Toilet Seat, 2
770-618	AQS 3IN1 RTS 4", WHITE, 1EA,2CV,Ret

SPECIFICATIONS

Description

Weight Capacity	300 lb(136 Kg)
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]