

## AquaSense Toilet Seat Riser with Hinge

HCPCS:



FIGURE A



Designed to help getting on and off the toilet for people with bending or sitting difficulties. Facilitates cleaning and hygiene for all users. Apart from providing independence, dignity and confidence while helping the user get on and off the toilet with minimal straining it will facilitate cleaning and hygiene for all users. When placed in the raised position, it will provide a non-obstructed area, for men.

- Hinge for ease of cleaning, used under the existing seat discreet, can go unnoticed in the users bathroom décor
- Adding 3 1/2" inches of height to your toilet making it easier and safer to sit or stand
- Made from PE that is both lightweight and durable to support up to 300 lbs
- Fits most standard toilets
- Feel good knowing that plastics parts of this product are made using recycled materials

### PRODUCTS

#### Description

770-630	AqS Hinged 3.5" Riser, White, 1EA,3CV,Re
770-631	AqS E.Hinge 3.5"Riser, White, 1EA,3CV,Re

### SPECIFICATIONS

#### Description

Weight Capacity	300 lb
-----------------	--------