

## Premium Plastic, Raised, Elongated Toilet Seat with Lock

HCPCS: E0244



FIGURE A



FIGURE B

- Designed for individuals who have difficulty sitting down or getting up from the toilet
- Adds 5" to height of toilet seat
- Tool-free arms are 21" apart and can be removed or added as needed
- Wide opening in front and back for personal hygiene
- Heavy-duty molded plastic construction provides additional strength and durability
- Standard locking mechanism ensures safety
- Lightweight and portable
- No tools required for installation
- Easy-to-clean
- Fits most elongated toilets
- Made using recycled materials

### PRODUCTS

#### Description

12013	Raised Toilet Seat w/Arms,1/cs, RTL
-------	-------------------------------------

### SPECIFICATIONS

#### Description

Warranty	Limited Lifetime
Seat Dimensions	17" (D) x 15.5" (W) x 5" (H)
Weight Capacity	300 lb
Carton Shipping Weight	6.5 lb
Weight	4.5 lb
Width Between Arms	21"