

Bariatric Transfer Bench

HCPCS: E0248



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Tool-free assembly of back, legs and arm
- A-frame provides additional stability
- Height adjusts in ½" increments with unique "Dual Column" extension legs
- All transfer benches are reversible to accommodate any bathroom
- Euro-style clip ensures a secure fit (Figure A)
- No exposed hardware to injure patient
- Large, locking suction cups for added safety (Figure B)
- Weight Capacity 500 lbs
- Made using recycled materials

PRODUCTS

Description

12025KD-1	Transfer Bench, Heavy Duty, KD, 1/cs
-----------	--------------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Seat Dimensions	26.5" (W) x 16.5" (D)
Weight Capacity	500 lb
Carton Shipping Weight	27 lb
Weight	14 lb