

Splash Defense Transfer Bench with Curtain Guard Protection

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G




FIGURE H

- New "U" shape cut-out design enables the shower curtain liner to stay within the bath tub and avoiding water spillage outside onto the bathroom floor
- Reversible to accommodate any bathroom
- Tool-free assembly of back, legs and arm
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs
- 400 lb weight capacity
- A-frame construction provides stability
- Durable blow-molded plastic bench and backrest
- Pinch-free lever allows for push pins to be depressed without pinching fingers
- Made using recycled materials

PRODUCTS

Description

RTL12032KDR	Splash Defense Transfer Bench
-------------	-------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Seat To Floor Height	17.5 22.5
Dimensions	Fits Within Tubs: Maximum Height 15.5, Minimum Bottom Width 18.5, Minimum Top Width 23, Maximum Edge Width 5.5
Height Adjustment	1/2" increments