

P903 Hand Pendant

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Controlled access to advanced functions using magnetic fob swipe to lock/unlock High/Low and Trendelenburg features.
- Easy control for patients with limited mobility thanks to a lightweight pendant and intuitive button layout.
- Convenient customization of bed positions with full control over head, foot, and overall height adjustments.
- Comfortable nighttime use enabled by backlit buttons for clear operation in low-light conditions.
- Enhanced safety through a lock feature that prevents accidental or unintended adjustments.
- Reliable long-term performance with a durable, water-resistant casing that withstands spills and drops.
- Quick access provided by an integrated clip for easy bedside mounting.
- Flexible positioning from any location with an extra-long, tangle-free cord.

PRODUCTS

Description

P903HP	P903 Hand Pendant, MagLockout, 1/ea
--------	-------------------------------------