

## Bariatric Aluminum Rehab Shower Commode Chair

HCPCS:



- Versatile 3-in-1 design enables this model to easily be positioned over a standard toilet, used as a portable self-contained commode or shower chair.
- Simply and quickly adjust the height to achieve the ideal seated position.
- Padded flip-back armrests provide superior support and create a comfortable, non-slip surface for a safe and secure grip when sitting or standing.
- Easy to empty commode bucket to reduce splashes and spills. Use commode liners for even quicker clean up and less mess!
- 5" locking casters allow for extra security when the unit is not in motion.
- Durable and easy-to-clean padded seat (22" x 17") and backrest.
- Seat opening allows for unobstructed access to the perineal area for easy and effective cleansing.
- Swing-away detachable footrests with height and angle adjustments.
- Both lightweight and durable the rust-resistant aluminum frame supports up to 500 lbs.
- Simple, tool-free assembly.

### PRODUCTS

### Description

NRS185008	Bari Shower&Commode Chair,1/ea
-----------	--------------------------------

### SPECIFICATIONS

### Description

Warranty	Limited Lifetime
Seat Dimensions	22" (W) x 17.1" (H)
Weight Capacity	500 lb
Seat To Floor Height	20.5" - 22.4"
Overall Width	26"