

Bariatric Folding Commode

HCPCS: E0168



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Tool-free assembly in 3 easy steps
- Width Between Legs: 29.25"
- Easily opens and folds for storage and transportation
- Durable plastic 16 ½" deep snap-on seat with lid
- Plastic armrests for added comfort
- Removable tool-free back
- Comes complete with 12 qt commode bucket with carry handle, cover and splash shield
- Weight Capacity 650 lbs
- Made using recycled materials

PRODUCTS

Description

| | |
|----------|--|
| 11117N-1 | FoldingBariatric Commode (DeepSeat),1/cs |
|----------|--|

SPECIFICATIONS

Description

| | |
|------------------------|------------------------|
| Warranty | Limited Lifetime |
| Seat Dimensions | 13.75" (W) x 16.5" (D) |
| Weight Capacity | 650 lb |
| Height | 15.5" - 22" |
| Carton Shipping Weight | 25 lb |
| Weight | 19 lb |