

Deluxe Steel Drop Arm Commode

HCPCS: E0165,E0165



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E



FIGURE F



FIGURE G



FIGURE H

Drop arm for ease of transfers, lightweight commode with increased width between arms.

PRODUCTS

PRODUCTS	Description
11125KD-1	Commode Drop Arm,Knock Down, 1/cs

SPECIFICATIONS

SPECIFICATIONS	Description
Product Weight	18 lb
Seat Dimensions	13.75" - 16.5"
Weight Capacity	300 lb
HCPCS	E0165
Height Adjustment	20" - 24"
Width Between Arms	19.75"
CA Prop 65	no