

Folding Steel Commode



HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Easily opens and folds
- Folds flat for convenient storage and transportation
- Comes complete with 7.5 qt commode bucket with carry handle, cover and splash shield
- Durable plastic snap-on seat and lid
- Assembles without pushing buttons
- Grey powder-coated steel-welded construction increases strength and durability
- Plastic arms for added comfort
- Made using recycled materials

PRODUCTS

Description

11148N-4	Commode,Folding,with DeepSeat,4/cs
----------	------------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Seat Dimensions	13.25" (W) x 16.5" (D)
Weight Capacity	350 lb
Height	15.5" - 21.75"
Carton Shipping Weight	59 lb
Weight	13 lb