

Tub Stand

HCPCS: E1399-1,E1399-1



- Slip-resistant legs
- Legs fold flat for storage
- Quick-release locks allow chair to be easily attached or removed
- Height adjustable 5" to 10"

PRODUCTS	Description
OT 8010	Otter Tub Stand, 1/cs

SPECIFICATIONS	Description
Weight Capacity	200 lbs
Carton Shipping Weight	11 lbs
Weight	8 lbs

PARTS	Description
OT TIPS	Replacement Tips,Otter, 4 ea/BX, BX
OT 8022	Caster Frame Insert, OT 8020, 1/ea, EA
OT 8021	Locking Caster for OT 8020,1/ea, EA
OT 8023	EndCapOtterShowerStand,1ea, EA
OT 8024	Replacement Screws for OT 8020, 4/bag, BAG
OT 0022	Release Lever, OT8010/OT8020, 1/ea, EA
OT 0023	Release Lever Assembly, OT 8020,1/ea, EA