

## Bariatric Aluminum Folding Cane, Height Adjustable

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Cane folds into four convenient parts for an easy storage
- Attractive wood handle with brass collar
- Manufactured with sturdy, extruded aluminum tubing
- Comes with plastic clip to hang cane when folded
- Handle height adjusts from 33" to 37"
- Weight Capacity 500 lb

### PRODUCTS

### Description

RTL10304HD	Cane,HD,Folding,Black,1EA,12/CS
------------	---------------------------------

### SPECIFICATIONS

### Description

Warranty	Limited Lifetime
Handle Height	33" - 37"
Weight Capacity	500 lb