

Nitro Aluminum Rollator, 10" Casters

HCPCS:




FIGURE A



FIGURE B




FIGURE C




FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

• Travel in style with the sleek design of the Euro-style Nitro Rollator. This rolling walker is cleverly engineered with a comfortable seat, convenient storage bag, adjustable height handles, and large 10 front casters for seamless rolling. With the brake cables inside the lightweight frame, the rollator walker provides a safe experience when you're out and about. With the ability to fold side-by-side, the Nitro Rollator is easily transported and stored.

- Adjustable contoured backrest
- Cross-brace design for added safety
- Clip hold and side-by-side folding
- Large and convenient storage bag included
- Weight capacity: 300 lb
- Weighs only 18 lb., easy-to-transport

PRODUCTS

Description

RTL10266	Nitro Rollator, Red Frame, 1/ea
RTL10266BK	Nitro Rollator, Black frame,1EA

SPECIFICATIONS

Description

Seat Dimensions	18" (W) x 10" (D)
Handle Height	33.5" - 38.25"
Weight Capacity	300 lb
Seat To Floor Height	20.5"
5 Year Limited Warranty	Brake and Handle Housing
Carton Shipping Weight	21.5 lb
Weight	18 lb

ACCESSORIES

Description

10266-CH	Nitro Cup Holder, 1/ea, EA
10266-BRAKE ENH	Nitro Slow Down Brake, ENH, 1PR, 25/CS, EA